

## RYCOLAB



The Film Free Shrink Tester determines the degree of unrestrained linear thermal shrinkage at given temperatures of plastic film and sheeting in fluid medium.

### Technical specifications

- Specimen Size  $\leq 160$  mm x 160 mm
- Temperature Range Room temperature ~ 200°C
- Accuracy  $\pm 0.3^\circ\text{C}$
- Power Supply AC 220V 50Hz
- Dimension 370 \* 440 \* 280 (LxWxH mm)

### Features

- P.I.D. temperature control
- High precision timer
- Alarming function upon test finish
- Stable heating of fluid medium
- Standard square metal mesh and free shrink holder

### Applicable standards

ASTM D2732