

Electrical circular sample cutter

Product code [RL-ECC-A](#)

RYCOLAB

Usage

For preparing 100 cm² circular samples on paper, cardboard, textile and other materials. For grammage or other properties determination.

Also available in smaller sizes: 50cm² or less.
Contact us for more information about the smaller versions.

Specifications

- Electric drive.
- Cutting thickness up to 7 mm.
- Autonomous and portable use. Possibility of working with batteries or electric supply feeding.
- Standard cutting area 100cm².
- Cutting by means of three blades simultaneously.
- Fast, clean and safe cut.
- CE mark.

Included accessories

- Power source.
- Allen wrench to replace the blades.
- Spare blades set.
- Support mat for cutting.

Connections

Electric : Rechargeable battery

Dimensions and weight

Only equipment: 210 x 210 x 160 (W x L x H) / 8 kg

With packaging box: 300 x 300 x 250 (W x L x H) / 10kg

