


## RYCOLAB



### Usage

For the preparation of several samples simultaneously.

### Device description

This apparatus allows to cut several strips simultaneously (edges being parallel). The precision is lower than 0.1 mm and the length depends on the sample. The samples are cut out by small wheels. There is no risk of fingers cutting.


### Models

Exist in various dimensions

Variable length

X= Number of strips

Y= Width of strips



Y	X	
1/2	1 to 5	12.7 mm width
15	1 to 5	15 mm width
25	1 to 5	25 mm width
01	1 to 5	25.4 mm width
30	1 to 5	30 mm width
50	1 to 4	50 mm width
100	1 to 2	100 mm width