

Disintegrator

Product code [RL-DSG-A](#)

RYCOLAB



Usage

For disintegration of pulp suspensions acc. to different standards.
Lightweight pot made from acrylic glass or heavy stainless steel pot

Applicable standards

ISO 5263-1 / SCAN C18 M2 / TAPPI T205 / PAPTAC C.6

Specification

- Made from robust stainless steel and aluminium
- Hinged head for easy take-out of the vessel
- Safety interlock
- Motor: 370 W
- Digital display of the revolutions
- Agitator speed according to standard: 2,980 +- 30 rpm applicable for TAPPI and
- ISO

Delivery content

- Disintegrator
- Acryl pot
- Measuring gauge
- Calibration certificate
- Operation manual
- Connection cable

Connections

Electric: 230 V- 50

Weight and dimensions

270 x 400 x 460 mm (WxLxH) / 36 Kg (Only Machine)

350 x 500 x 600 mm (WxLxH) / 56 Kg (With packaging box)